Metabolite Profiling Analysis of Conventional and Non-Conventional Extraction Methods on Secondary Metabolite from Peperomia pellucida (L.) Kunth using UPLC-QToF-MS/MS System

by Islamudin Ahmad

Submission date: 02-Sep-2019 08:01AM (UTC-0700) Submission ID: 1166260038 File name: JYoungPharm_10_2_s40.pdf (660.55K) Word count: 5006 Character count: 26493 www.jyoungpharm.org | www.phcog.net

Metabolite Profiling Analysis of Conventional and Non-Conventional Extraction Methods on Secondary Metabolite from *Peperomia pellucida* (L.) Kunth using UPLC-QToF-MS/MS System

Islamudin Ahmad^{1,3}, Kamarza Mulia², Arry Yanuar³, Abdul Mun'im^{3*}

¹Laboratory of Pharmaceutical Research and Development of FARMAKA TROPIS, Faculty of Pharmacy, Mulawarman University, Samarinda, East Kalimantan, 75119, INDONESIA.

²Department of Chemical Engineering, Faculty of Engineering, Universitas Indonesia, Depok, West Java, 16424, INDONESIA. ³Department of Pharmaceutical Sciences, Faculty of Pharmacy, Universitas Indonesia, Depok, West Java, 16424, INDONESIA.

ABSTRACT

Objective: The purpose of this study was to observe the difference of extraction method (both conventional and non-conventional) based on metabolite profile using UPLC-QToF-MS/MS system from Peperomia pellucida (L.) Kunth. Methods: Dried samples were extracted using the conventional maceration method and the optimum ofionic liquid-based microwave-assisted extraction (IL-MAE) methods.Metabolite profiling was performed using UPLC-QToF-MS/MS system with some modifications adjusted to the instrument condition. The data was analyzed using Masslynx 4.1 software. Results: based on the results, there were differences on metabolite profiling from both conventional and non-conventional extraction methods that was extraction method using organic solvent and ionic liquid solvent ([BMIM]BF, at the optimum condition. The extract obtained using IL-MAE method had a peak depth with a well-separated Rt (retention time) value ranging from 0.5 to 7.5 min which means that the extracted compound was from polar to nonpolar properties. The extract obtained using maceration method, the peak spread on a separate Rt ranges from 2.5 to 7.5 min. Also, both extracts obtained have different area under curve (AUC) values

with amount total of 22285 (ILMAE) and 12679 (maceration), respectively, and showed the ILMAE twice as large as that of the maceration. **Conclusion**: based on the results, conventional and non-conventional extraction method showed differences in metabolite profiling based on Rt value and mass spectrum m/z of each peak.

Key words: Peperomia pellucida (L) Kunth, Hyphenated system, Ionic liquid-based microwave-assisted extraction, Maceration, Metabolite profiling.

Correspondence

Abdul Mun'im, Department of Phar Maceutical Sciences, Faculty of Pharmacy, Universitas Indonesia, Depok, West Java, 16424, INDONESIA.

Phone: +62-852-1610-4550

Email: munimabdoel@gmail.com DOI: 10.5530/jyp.2018.2s.8

INTRODUCTION

Peperomia pellucida (L.) Kunth herbs (Piperaceae family) has been known to contain secondary metabolites such as terpenoids, alkaloids, saponins, dan polyphenols.¹ Some of them have been isolated (particularly polyphenolic groups) including pellucidin A,² peperomins, secolignan,³ xanthone patuloside A,⁴ dillapiole,⁵ chromene,⁶ and quercetin.⁷ Also, the herbs have also been reported to have potential pharmacological properties.⁸

The bioactive compound from natural products mainly from plants source can be extracted by conventional (including maceration, reflux, soxhlet, percolation)⁹ and non-conventional extraction methods (such as microwave-assisted extraction, ultrasonic-assisted extraction, supercritical fluid extraction, and so on).¹⁰ However, the use of extraction methods is strongly influenced by the type of solvent used.

lonic liquid based microwave assisted extraction (IL-MAE) continued to be developed and applied in the extraction secondary metabolites from natural products. This method was developed based on the selection of solvents with a green chemistry principle approach in exploring the metabolite constituent of natural products. The ionic liquid is chosen as solvent because it has the flexibility and potential to replace conventional organic solvents that are volatile, flammable, and toxic.¹¹ Application of ionic liquid as a solvent to extracting secondary metabolite and optimization of the IL-MAE method of polyphenolics content extraction from this herb has been performed in previous studies.¹²⁻¹³

Metabolite profiling is intended to distinguish between the known compound (dereplication) and new molecule directly present in the crude extract.¹⁴ Metabolite profiling in crude extracts is not an easy work as natural materials show a very diverse form of structure. In some compounds, the atomic sequence and its stereochemical orientation must be elucidated by de novo. Consequently, the ability of a single analytical technique in profiling metabolite content in the crude extract to date does not exist.¹⁵ Innovative development strategies from metabolite profiling in crude extracts where the researcher can be facilitated by the hyphenated system found during the last decade. The potential for chemical screening strategies has been significantly enhanced by the recent development of the hyphenated technique, which is capable of generating efficient metabolite separation, along with valuable structural information both online and at-line.¹⁶

One of the instruments with the hyphenated system used to analyze metabolite profiling with high accuracy and validity¹⁷⁻¹⁸ from crude extract is UPLC-QToF-MS/MS tandem system.¹⁹ This method has been

This is an open access article distributed under the terms of the Creative Commons Attribution-NonCommercial-ShareAlike 4.0 License, which allows others to remix, tweak, and build upon the work non-commercially, as long as the author is credited and the new creations are licensed under the identical terms.

Journal of Young Pharmacists, Vol 10, Issue 2 (Suppl), Apr-Jun, 2018

Original Article

successfully used for metabolite profiling analysis in some medicinal plants such as metabolite profiling of bioactive compounds from *Bidens pilosa*,¹⁴ metabolites profiling secondary metabolites of *Arabidopsis arabika*,²⁰ Analysis of steviol and its glycoside from *Stevia rebaudiana* leaves as a commercial sweetener,²¹ metabolite profiling of polyphenols from *Vaccinium berries*,²² While the metabolite profiling analysis of *P. pellucida* herb (using both conventional and non-conventional extraction) has not been reported.

Active compound isolation from *P. pelucida* herb requires specific techniques to obtain the compound optimally. Development of extraction methods (both conventional and non-conventional) and metabolite profiling with hyphenated systems to achieve the target compounds. This study aimed to know the difference of extraction method (both conventional maceration and non-conventional IL-MAE) based on metabolite profile analysis using UPLC-QToF-MS/MS system.

MATERIALS AND METHODS

Materials

Wild herb materials of *P. pellucida* were collected from the oil palm plantation at Baras sub-district, North Mamuju, West Sulawesi, Indonesia, and were identified at the Herbarium Bogoriense, Bogor, West Java, Indonesia. The fresh sample was washed, and dried at 50–60°C using drying oven, and the dried sample was powdered using a grinder. The powder obtained was stored at a cold temperature until analysis.

Material and General Equipments

The materials used in this study include n-hexane, ethyl acetate, aquadest (PT. Smart Lab Indonesia, Indonesia), 1-butyl-3-methylimidazolium tetrafluoroborate (Cheng Jie Chemical Co. LTD, China). Acetonitrile (LiChrosolv* for HPLC), Ethanol (Chromosolv* for HPLC), formic acid (Merck, Germany). The general equipment used in this study were Modena Microwave 900 Watt (Buono-MV3002), Rotary Evaporator (Buchi, Germany), UPLC-QTOF-MS/MS System (Waters, USA), as well as equipment for maceration.

Extraction procedure

The extraction process was conducted using conventional (maceration method) and non-conventional (ILMAE method) as follows:

1. Conventional Maceration Method

Dried sample (1 kg) was macerated using n-hexane solvent for 24 h, and extraction process was performed for 4-6 times. Next, the residual of samples was re-macerated using ethyl acetate with the same condition. The extract solution was evaporated using a rotary evaporator to obtain a dried extract.¹²

2. Non-conventional IL-MAE Method

For extraction process with the optimum ILMAE method was performed based on the previous study.^{12,23} including microwave power, extraction time, ionic liquid concentration, and liquid-solid ratio had been involved. Response surface methodology and Box Behnken design were used to obtain predictive model (multivariate quadratic regression equation Briefly,30 grams of dried sample was extracted using the optimum condition of IL-MAE method with [BMIM]BF₄ as a solvent. The optimum condition of IL-MAE method was performed using ionic liquid concentration of 0.79 mol/l, the liquid-solid ratio of 11 ml/g sample, extraction time of 20 min, and microwave power of 30 %Watt.²³ For analysis, the extraction process was performed four times.

Metabolite profiling analysis

For determination of metabolite profiling was performed using UPLC-QToF-MS/MS System according to the literature,^{14,17,22,24} with

Journal of Young Pharmacists, Vol 10, Issue 2 (Suppl), Apr-Jun, 2018

some modification adjusting to the instrument condition. Briefly, liquid chromatographic (LC) separation was performed using Xevo G2 ACQUITY with BEH C₁₈ column 1.7 µm, 2.1 x 50 mm (Waters Corp, Milford, MA, USA). Autosampler Tray and columns temperature was maintained at 40°C, the mobile phase consisting of water containing 0.1% formic acid (A) and acetonitrile containing 0.1 formic acids (B) was used at a flow rate 0.3 ml/min. The elution gradient can be seen in Table 1. Data acquisition was processed with Masslynx 4.1 software. Extract sample (10 mg) was weighed and dissolved in 100 ml methanol (LC-MS grade). One ml of the sample solution was taken into the Effendorf tube and was centrifuged at 13000 rpm for 10 min. The supernatant was separated into another Effendorf tube. The sample injection volume was set at 5 µl and the total running time was 10 min. The UPLC output system was connected to the mass spectrometer via the ESI interface and was operated in the positive ion detection mode. Nitrogen was used for desolvation of 500 L/h; the cone gas was set at 16 L/h. The capillary voltage was set at 3 kV, source temperature 110°C, and desolvation temperature 300°C.

RESULTS

For the results test of profiling metabolite using the same gradient conditions that were the mobile phase with the ratio between aquadest and acetonitrile (containing 0.1% formic acid) from 95:5 to 5:95 for 9 min as shown in Table 1.

As can be seen in Figure 1 and 2 demonstrated Total Ion Chromatogram (TIC) with peak smoothing and peak area integration (A) and TIC with peak smoothing (B) that showed difference of metabolite profiling from both method that was extraction method using organic solvent and ionic liquid solvent ([BMIM]BF₄) at the optimum condition of the instrument. In Figure 1 a peak depth with a well-separated Rt (retention time) value ranging from 0.5 to 7.5 min was shown, which means that the extracted compound was from polar to nonpolar properties. In Figure 2, the peak spread on a separate Rt ranges from 2.5 to 7.5 min. While in the min of 0.5 to 2.4 there was no peak. It indicates that the compound extracted by maceration method using ethyl acetate have nonpolar properties.

From both extracts based on different extraction methods (i.e., maceration and IL-MAE) the different result of metabolite profiling was obtained. Shown in Figure 1 which was the result of metabolite profiling of the extract obtained by IL-MAE, there was a peak especially in 0.4 to 3.0 min which was not present in the extract obtained by maceration (as can be seen in Figure 2). Whereas in Figure 2 there were peaks, especially at strong peak with Rt of 4.52 – 5.60 min while peak on metabolite profiling results from the extract obtained by IL-MAE show peak which was weak even not exist. While in the peak with Rt 6.20 there were on both extracts. Also, both extracts obtained has different AUC values in which extracts

Table 1: The elution gradient for metabolite profiling analysis using UPLC- QToF-MS/MS System (Waters).

				-
No.	Time (Minute)	%A	%B	
1	0	95	5	l
2	1	95	5	
3	6	0	100	
4	7	0	100	
5	7.5	95	5	
6	9	95	5	
				1

Where A is the mobile phase consisting of water containing 0.1% formic acid and B is the mobile phase consisting of acetonitrile containing 0.1% formic acid.





Figure 1: Peaks data of metabolite profiling analysis from extracts were obtained by IL-MAE method (non-conventional method) with [BMIM] BF_aas a solvent using UPLC- QToF-MS/MS System. Where A is Total ion chromatogram (TIC) with peak smoothing and peak area

integration; B is Total ion chromatogram with smoothing.

Figure 2: Peaks data of metabolite profiling analysis from extracts were obtained by maceration method (conventional method) with ethyl acetate as a solvent using UPLC- QToF-MS/MS System. Where A is Total ion chromatogram (TIC) with peak smoothing and peak area integration; B is Total ion chromatogram with smoothing.

Table 2: Results of UPLC-QToF-MS/MS (m/z value) spectrum of extract was obtained using IL-MAE method (non-conventional method).

No	Retention Time (Rt)	min	m/zvalue		
1	0.4 - 2.3	0.468	139.1816; 140.1761; 141.1166; 323.3341; 507.4559; 691.6055; 875.7525; 1059.8871; 1244.0045; 1428.1177		
		0.625	139.1782; 140.1761; 364.3413; 365.3409; 366.3379; 591.4858; 626.9603; 817.6363; 922.6903; 1043.7958; 1268.9130; 1496.0017; 1496.0001; 1496.0001; 1496.0001; 1496.0001; 1496.0001; 1496.0000; 1496.00		
		0.730	$139.1768; 140.1758; 365.3379; 410.1645; 591.4912; 626.9582; 817.6421; 922.6833; 1043.7909; 1268.9534; 1495.0270 \\ 1268.954; 1495.0270 \\ 1268.9545; 1495.0270 \\ 1268.9545; 1495.0270 \\ 1268.9545; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0270 \\ 1268.9555; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.0255; 1495.025; 1495.0$		
		2.082	211.2456; 212.2475; 251.2299; 354.1744; 486.2294		
		2.232	211.2516; 212.2471; 251.237; 355.1892; 565.2423; 647.2996; 818.2544; 1054.0790		
2	2.3 - 2.9	2.436	223.2466; 245.2317; 252.2227; 340.1905; 565.2446; 603.3232		
		2.490	171.1478; 245.2336; 246.2336; 251.2426; 433.2215; 481.2777; 595.3393		
		2.669	195.1249; 196.1307; 283.2699; 447.2441; 593.3315; 685.3773		
		2.715	195.1291; 196.1336; 402.2793; 477.2263; 522.3278; 687.3209		
		2.805	179.1615; 219.1680; 308.3018; 403.2110; 566.4335; 724.3237; 769.3845; 911.3472; 1027.4069; 1174.3456; 1494.6885; 1494.685;		
3	2.9 - 3.5	2.919	201.2108; 219.2210; 235.2005; 323.3064; 407.2352; 430.2371; 531.3644; 650.2902; 847.4064; 899.3586		
		3.091	151.1069; 209.1243; 211.1613; 276.1473; 385.2653; 407.2422; 515.3213; 647.4405; 733.3910		
		3.109	177.1185; 207.1454; 314.2209; 331.2057; 407.2459; 548.3253; 549.3280; 671.3453; 783.3347; 883.4658		
		3.348	207.1514; 225.1598; 305.2462; 421.2552; 529.3532; 549.3170; 685.3718; 685.3718; 773.4330; 957.6433; 1175.5614; 568.3718; 568		
		3.431	195.1424; 271.1469; 313.1953; 373.2277; 375.2440; 478.3654; 593.3487		
4	3.5 - 4.5	3.685	212.2958; 304.3846; 313.1908; 318.3595; 529.3406; 625.3709; 754.4578; 841.4797; 975.5499		
		3.756	$207.1587; 235.1687; 236.1685; 403.2746; 404.2841; 529.3494; 753.4700; 827.4679; 952.3814; 1050.5581; 1189.4120; 1438.4369 \\ 1060000000000000000000000000000000000$		
		3.900	225.2053; 384.4238; 285.4135; 332.4052; 525.2990; 561.3804; 739.4343; 817.4390; 949.5057; 1075.5227		
		3.953	$300.1483; 315.1792; 316.1764; 332.3993; 466.3031; 629.3246; 739.4461; 785.4274; 963.5539; 999.5621; 1111.6162 \\ 100000000000000000000000000000000000$		
		4.129	195.1401; 312.4749; 313.4536; 360.4524; 581.3390; 653.3799; 739.4486; 869.4777; 986.4565; 1124.5309		
		4.218	195.1435; 322.4378; 329.1913; 330.1898; 331.2008; 441.4445; 619.3845; 739.4498; 949.5218; 1174.6296; 1286.6456; 1286.6466; 1286.6456; 1		
5		4.959	$149.0787; 277.2780; 301.2099; 522.4457; 623.3851; 766.6262; 898.7619; 987.7861; 1174.6445; 1313.7004; 1495.8024 \\ 149.0787; 277.2780; 301.2099; 522.4457; 623.3851; 766.6262; 898.7619; 987.7861; 1174.6445; 1313.7004; 1495.8024 \\ 149.0787; 277.2780; 301.2099; 522.4457; 623.3851; 766.6262; 898.7619; 987.7861; 1174.6445; 1313.7004; 1495.8024 \\ 149.0787; 277.2780; 301.2099; 522.4457; 623.3851; 766.6262; 898.7619; 987.7861; 1174.6445; 1313.7004; 1495.8024 \\ 149.0787; 277.2780; 301.2099; 522.4457; 623.3851; 766.6262; 898.7619; 987.7861; 1174.6445; 1313.7004; 1495.8024 \\ 149.0787; 149.0787; 149.0864; 149.0866; 149.0866; 149.0866; 149.0866; 149.0866; 149.0866; 140.0866; 140.0866; 140.0866; 140.0866; 140.0866; 140.0866; 140.0866; 140.0866; 140.08666; 140.0866; 140.0866; 140.0866; 140.0866; 140.0866; 140.0866; 140$		
		5.714	139.1722; 251.1402; 308.3557; 309.3383; 419.3854; 565.3870; 581.3549; 637.4206; 638.4244; 639.4299; 858.8097; 896.7033; 1053.4819; 1272.7568; 1460.0713		
		6.218	$149.0884; 150.0798; \ 391.3789; \ 413.3726; \ 429.3433; \ 457.3602; \ 638.7296; \ 803.6777; \ 872.8376; \ 1019.749$		
		7.155	139.1703; 149.0715; 370.3681; 413.3850; 607.5057; 663.5709; 843.6861; 887.7139; 888.7165; 889.7233; 890.7280; 985.9030; 1326.0479		
		7.206	139.1727; 370.3823; 413.3847; 677.5908; 843.6845; 871.7186; 872.7230; 873.7269; 874.7251; 1090.8916; 1155.9614; 1342.0386; 1470.1021		

Journal of Young Pharmacists, Vol 10, Issue 2 (Suppl), Apr-Jun, 2018

	No	Retention Time (Rt)	Min	m/zvalue
1	1	3.34 - 4.3	3.635	181.1697; 219.1656; 304.3615; 387.3539; 753.4727; 814.4013; 924.4348; 1206.4136; 1285.6440
			3.756	219.1701; 235.11686; 236.1662; 329.2194; 403.2748; 427.2461; 632.3526; 827.4625; 81.4370; 976.5239
			3.950	219.1695; 315.1745; 316.1740; 427.2531; 553.3780; 682.2105; 773.4325
			4.140	219.1689; 312.4491; 313.4427; 315.1689; 427.2712; 643.3893; 720.4321; 947.6596; 983.5687; 1143.5848; 1324.7644; 1436.8271
			4.272	219.1684; 221.1547; 329.1894; 330.1904; 427.2780; 627.3904; 711.3359; 803.5481; 967.4610
	2	4.3 - 5.3	4.400	219.1678; 221.1811; 299.1703; 304.2689; 411.2731; 627.3802; 699.4662; 737.5193; 880.5945
			4.512	221.1875; 222.1843; 411.2766; 553.4916; 627.3677; 795.5170; 833.5162; 981.5439; 1091.9723
			4.941	203.2382; 219.1720; 221.1771; 387.2785; 419.3132; 639.3691; 641.3979; 810.6937; 939.6212; 1031.7057
			5.213	219.1744; 275.2413; 387.2763; 607.3751; 610.3864; 611.3885; 669.4666; 879.4844; 1022.9031; 1217.7123; 1269.7792; 1494.8412
			5.277	219.1739; 386.3646; 607.3882; 609.4003; 635.3800; 626.3835; 627.3950; 782.6851; 829.6377; 1029.4808; 1233.7205; 1357.7954; 1434.646
	3	5.3 - 5.85	5.363	219.1739; 387.2882; 509.3809; 535.3666; 610.3857; 611.3904; 782.6791; 879.5152; 1086.8999; 1218.17064; 1318.7473 + 100000000000000000000000000000000000
			5.438	219.1713; 412.3825; 535.3679; 593.3825; 594.3857; 609.3829; 625.4382; 782.6993; 1127.7306; 1185.7092; 1354.2301
			5.560	219.1709; 221.1738; 341.2921; 414.4041; 415.4006; 535.3690; 593.3873; 594.3924; 743.4485; 829.6261; 919.6580; 1197.0258
			5.714	203.2279; 219.1732; 308.3825; 309.3676; 535.3777; 593.3900; 637.4399; 829.6188; 919.6422; 1256.0709; 1313.7601
			5.785	219.1567; 334.4149; 335.3990; 356.3844; 607.4096; 689.6401; 813.6051; 903.6359
	4	5.85 - 7.5	5.907	203.2272; 219.1738; 220.1809; 348.4148; 349.4074; 387.2815; 611.5715; 612.5655; 829.6135; 905.6819; 1097.9205; 1188.9225; 1376.1067
			6.014	219.1734; 221.1704; 336.5163; 362.4325; 427.3952; 551.4983; 815.6162; 887.6426; 937.7065; 1042.7949; 1167.8545
			6.204	149.0832; 219.1691; 387.2775; 413.3725; 429.3481; 613.5667; 797.6502; 865.6614; 961.7382; 1065.7802
			6.491	219.1746; 335.3382; 387.2788; 429.4297; 613.5734; 797.6614; 905.7272; 921.7256; 922.7267; 923.7308; 1033.8214; 1159.9659
			7.141	219.1692; 370.3982; 541.3345; 609.3832; 871.7190; 903.7157; 904.7144; 905.7189; 906.7194; 1081.8163; 1275.9602; 1497.1693
			7.238	219.1591; 370.3939; 593.3926; 705.6019; 813.7011; 871.7213; 872.7220; 873.7233; 874.7306; 1093.7841; 1166.9430; 1296.0737; 1433.9531
			7.395	219.1620; 370.3899; 565.3514; 705.6183; 813.7059; 814.7114; 815.7146; 888.7054; 1093.8440; 1166.9338; 1412.1720

Table 3: Results of LC-MS (m/z value) spectrum of extract was obtained using maceration method (conventionalextraction method)." replaced with "Results of UPLC-QToF-MS/MS (m/z value) spectrum of extract was obtained using maceration method (conventional extraction method).

were extracted by IL-MAE and maceration methods with amount total of 22285 and 12679, respectively, where the extracts of IL- MAE was twice as large as that of the extracts of maceration method. replaced with "method. The m/z values of each peak based on different Rt values demonstrated in Table 2 for IL-MAE extract and Table 3 for ethyl acetate extract.

DISCUSSION

In this study, non-conventional extraction was performed using ILMAE method based on a combination of optimum extraction condition parameters, as well as 18.6 min extraction times, 30 %Watt microwave power, 11 ml/g liquid-solid ratio, and 0.79 ml/L [BMIM]BF₄ in water. [BMIM]BF₄ was chosen because based on the comparison of effectiveness in extracting the polyphenol class from this herb in the previous study, and the liquid ionic solvent was the best solvent to extract the target compound.²³ including microwave power, extraction time, ionic liquid concentration, and liquid-solid ratio had been involved. While the conventional extraction using the maceration method with n-hexane then followed by ethyl acetate. Each extract obtained was dried and stored in a cool place and prepared for analysis. Strategies for the analysis and identification of natural plant compounds with profiling metabolite approach using UPLC-QToF-MS/MS tandem system based on peaks at a particular retention time (Rt) as data and compared with the database,¹⁹ particularly natural product research. In metabolite profiling analysis, a very high resolution of chromatography was required to detect the metabolites of the associated matrix complex.²⁶ The reduction of particle size in chromatography allows higher separation efficiency and shorter analysis times.¹⁶

In this study, metabolite profiling analysis was performed only limited to the results of the investigation using UPLC-QTOF-MS/MS tandemsystem and obtained peak distribution from each extract at the same eluent condition based on total ion chromatogram (TIC) with peak smoothing and peak area integration. The m/z spectrum value in each peak indicates compound predictions according to the experimental molecular weight approach as can be seen in Table 2 and Table 3. This study does not use standards as controls, however, focusing only on the differences in the components present in each extract based on the Rt and m/z spectrum value.

Based on the results of metabolite profiling analysis, it was shown that the results of spectrum m/z gave different peaks on each extract. However, in this study, the m/z spectrum cannot be determined precisely given the absence of a comparison standard to be used as a benchmark on each

Journal of Young Pharmacists, Vol 10, Issue 2 (Suppl), Apr-Jun, 2018

peak produced. In addition, these profiles correlate with total polyphenolics content from extracts of each method.^{13,23} Therefore, only the optimum conditions obtained from the instrument used to separate the class of compounds contained in the extract based on the technique used, and the preliminary data related to the use of liquid ionic solvents as an alternative solution with a green chemical principles approach.

CONCLUSION

Based on the above results, metabolite profiling was determined and analyzed using UPLC- QToF-MS/MS System and both conventional and non-conventional extraction method showed differences in metabolite profiling based on Rt value and mass spectrum m/z of each peak.

ACKNOWLEDGEMENT

This study was funded by Faculty of Pharmacy, Mulawarman University via a grant "Hibah Fakultas Farmasi UNMUL 2017." Many thank Head Laboratory of Pharmacognosy-Phytochemistry, Faculty of Pharmacy, Universitas Indonesia for providing facilities.

CONFLICT OF INTEREST

The author declares that he has no conflict of interest.

ABBREVIATIONS

P. pellucida: *Peperomia pellucida* (L.) Kunth; **IL-MAE:** Ionic Liquid Microwave-Assisted Extraction; **[BMIM]BF4:** 1-butyl-3-methylimidazolium tetrafluoroborate; **LC:** liquid Chromatography; **HPLC:** High Performance Liquid Chromatography; **UPLC-QToF-MS/MS:** Ultra-Performance Liquid Chromatography quadrupole time of flight mas spectrometry; **Rt:** Retention time; **TIC:** Total Ion Chromatogram; **AUC:** Area Under Curve; **m/z:** mass per atomic charge.

SUMMARY

- There were differences on profiling metabolite from both conventional and non-conventional extraction methods that was extraction method using organic solvent and ionic liquid solvent ([BMIM]BF4) at the optimum condition.
- The extract obtained using IL-MAE method had a peak depth with a well-separated Rt (retention time) value ranging from 0.5 to 7.5 minutes which means that the extracted compound was from polar to non-polar properties.
- The extract obtained using maceration method, the peak spread on a separate Rt ranges from 2.5 to 7.5 minutes.
- Both extracts obtained have different area under curve (AUC) values with amount total of 22285 (IL-MAE) and 12679 (maceration), respectively, and showed the IL-MAE twice as large as that of the maceration.
- Conventional and non-conventional extraction method showed differences in metabolite profiling based on Rt value and mass spectrum m/z of each peak.

REFERENCES

- Mensah JK, Okoli RI, Turay AA, Ogie-Odia EA. Phytochemical analysis of medicinal plants used for the management of hypertension by Esan people of Edo State, Nigeria. Ethnobot Leafl. 2009;13(12):73-87.
- Verma RS, Padalia RC, Goswani P, Chauhan A. Essential oil composition of *Peperomia pellucida* (L.) from India. J Essent Oil Res. 2014;26(1):29-33. doi:10. 1080/10412905.2013.822431.
- Xu S, Li N, Ning MM, Zhou CH, Yang QR, Wang MW. Bioactive Compounds from Peperomia pellucida. J Nat Prod. 2006;69(2):247-50. doi:10.1021/np050457s.

- Khan A, Rahman M, Islam S. Isolation and Bioactivity of a Xanthone Glycoside from *Peperomia pellucida*. Life Sci Med Res. 2010;2010:1-10.
- Rojas-Martínez R, Arrieta J, Cruz-Antonio L, Arrieta-Baez D, Velázquez-Méndez AM, Sánchez-Mendoza ME. Dillapiole, isolated from *Peperomia pellucida*, shows gastroprotector activity against ethanol-induced gastric lesions in wistar rats. Molecules. 2013;18(9):11327-37. doi:10.3390/molecules180911327.
- Susilawati Y, Nugraha R, Muhtadi A, Soetardjo S, Supratman U. (S)-2-Methyl-2-(4-methylpent-3-enyl)-6-(propan-2-ylidene)-3,4,6,7-tetrahydropyrano[4,3-g] chromen-9(2H)-one. Molbank. 2015;2015(2):1-6. doi:10.3390/M855.
- Kurniawan A, Saputri FC, Ahmad I, Mun'im A. Isolation of angiotensin converting enzyme (ACE) inhibitory activity quercetin from *Peperomia pellucida*. Int J PharmTech Res. 2016;9(7):115-21.
- Mun'im A, Nurpriantia N, Setyaningsih R, Syahdi RR. Optimization of microwave-assisted extraction of active compounds, antioxidant activity, and angiotensin-converting enzyme (ACE) inhibitory activity from *Peperomia pellucida* (L.) Kunth. J Young Pharm. 2017;9(1):73-8. doi:10.5530/jyp.2017.
- 9. Hanani E. Analisis Fitokimia. Jakarta: Penerbit Buku Keokteran, EGC. 2015.
- Khoddami A, Wilkes M, Roberts T. Techniques for analysis of plant phenolic compounds. Molecules. 2013;18(3):2328-75. doi:10.3390/molecules18022328.
- Jessop PG, Jessop D, Fu D, Phan L. Solvatochromic parameters for solvents of interest in green chemistry. Green Chem. 2012;14(5):1245-59. doi:10.1039/ c2gc16670d.
- Ahmad I, Yanuar A, Mulia K, Mun'im A. Application of ionic liquid-based microwave-assisted extraction of the secondary metabolite from *Peperomia pellucida* (L) Kunth. Pharmacogn J. 2017;9(2):227-34. doi:10.5530/pj.2017.2.38.
- Ahmad I, Yanuar A, Mulia K, Mur'im A. Application of ionic liquid as a green solvent for polyphenolics content extraction of *Peperomia pellucida* (L) Kunth herb. J Young Pharm. 2017;9(4):S1-S4. doi:10.5530/jyp.2017.
- Chiang YM, Chuang DY, Wang SY, Kuo YH, Tsai PW, Shyur LF. Metabolite profiling and chemopreventive bioactivity of plant extracts from *Bidens pilosa*. J Ethnopharmacol. 2004;95(2-3):409-19. doi:10.1016/j.jep.2004.08.010.
- Kaufmann A. Combining UHPLC and high-resolution MS: A viable approach for the analysis of complex samples? TrAC - Trends Anal Chem. 2014;63:113-28. doi:10.1016/j.trac.2014.06.025.
- Chechinel-Filho V. Plants Bioactives and Drug Discovery: Principles, Practice, and Perspectives. Hoboken, New Jersey: John Wiley and Sons, Inc. 2012.
- Cuthbertson DJ, Johnson SR, Piljac-Žegarac J, Kappel J, Schäfer S, Wüst M, et al. Accurate mass-time tag library for LC/MS-based metabolite profiling of medicinal plants. Phytochemistry. 2013;91:187-97. doi:10.1016/j.phytochem.2013.02.018.
- González RR, Fernández RF, Luis J, Vidal JL, Frenich AG, Pérez ML. Development and validation of an ultra-high performance liquid chromatography-tandem mass-spectrometry (UHPLC-MS/MS) method for the simultaneous determination of neurotransmitters in rat brain samples. J Neurosci Methods J. 2011;198(2):187-94. doi:10.1016/j.jneumeth.2011.03.023.
- Wolfender JL, Marti G, Thomas A, Bertrand S. Current approaches and challenges for the metabolite profiling of complex natural extracts. J Chromatogr A. 2015;1382:136-64. doi:10.1016/j.chroma.2014.10.091.
- D'Auria JC, Gershenzon J. The secondary metabolism of Arabidopsis thaliana: growing like a weed. Curr Opin Plant Biol. 2005;8(3):308-16. doi:10.1016/j. pbi.2005.03.012.
- Gardana C, Scaglianti M, Simonetti P. Evaluation of steviol and its glycosides in Stevia rebaudiana leaves and commercial sweetener by ultra-high-performance liquid chromatography-mass spectrometry. J Chromatogr A. 2010;1217(9):1463-70. doi:10.1016/j.chroma.2009.12.036.
- Prencipe FP, Bruni R, Guerrini A, Rossi D, Benvenuti S, Pellati F. Metabolite profiling of polyphenols in *Vaccinium berries* and determination of their chemopreventive properties. J Pharm Biomed Anal. 2014;89:257-67. doi:10.1016/j.jpba.2013.11.016.
- Ahmad I, Yanuar A, Mulia K, Mun'im A. Optimization of ionic liquid-based microwave-assisted extraction of polyphenolic content from *Peperomia pellucida* (L) Kunth using response surface methodology. Asian Pac J Trop Biomed. 2017;7(7):660-5. doi:10.1016/j.apjtb.2017.06.010.
- Grata E, Boccard J, Guillarme D, et al. UPLC TOF-MS for plant metabolomics : A sequential approach for wound marker analysis in Arabidopsis thaliana & J Chromatogr BJ. 2008;871(2):261-70. doi:10.1016/j.jchromb.2008.04.021.
- Hafidz KA, Puspitasari N, Azminah, Yanuar A, Artha Y, Mun'im A. HMG-CoA reductase inhibitory activity of *Gnetum gnemon* seed extract and identification of potential inhibitors for lowering cholesterol level. J Young Pharm. 2017;9(4):559-65. doi:10.5530/jyp.2017.
- Creek DJ, Dunn WB, Fiehn O, *et al.* Metabolite identification: are you sure? And how do your peers gauge your confidence? Metabolomics. 2014;10(3):350-3. doi:10.1007/s11306-014-0656-8.

Article History: Submission Date : 29-11-2017 ; Revised Date : 28-01-2018; Acceptance Date : 14-02-2018. Cite this article: Ahmad I, Mulia K, Yanuar A, Mun'im A. Metabolite Profiling Analysis of Conventional and Non-Conventional Extraction Methods on Secondary Metabolite from *Peperomia pellucida* (L.) Kunth Using .UPLC-QToF-MS/MS System. J Young Pharm. 2018;10(2)Suppl:s40-s44. Metabolite Profiling Analysis of Conventional and Non-Conventional Extraction Methods on Secondary Metabolite from Peperomia pellucida (L.) Kunth using UPLC-QToF-MS/MS System

ORIGINALITY REPORT 4% 4% 6% PUBLICATIONS SIMILARITY INDEX INTERNET SOURCES STUDENT PAPERS **PRIMARY SOURCES** www.thieme-connect.com % Internet Source Marina Ramos, Alfonso Jiménez, María Carmen 1% 2 Garrigós. "II-based advanced techniques for the extraction of value-added compounds from natural sources and food by-products", TrAC Trends in Analytical Chemistry, 2019 Publication Submitted to Yeungnam University <1% 3 Student Paper <1% www.ncbi.nlm.nih.gov 4 Internet Source <1% Armando, Y.P.. "Bioequivalence assay between 5 orally disintegrating and conventional tablet formulations in healthy volunteers", International Journal of Pharmaceutics, 20090121 **Publication**

6	buleria.unileon.es Internet Source	<1%
7	Submitted to University of Leicester Student Paper	<1%
8	www.freepatentsonline.com	<1%
9	Wolfender, Jean-Luc, Guillaume Marti, Aurélien Thomas, and Samuel Bertrand. "Current approaches and challenges for the metabolite profiling of complex natural extracts", Journal of Chromatography A, 2015. Publication	<1%
10	westminsterresearch.westminster.ac.uk	<1%
11	ujcontent.uj.ac.za Internet Source	<1%
12	Winny Routray, Valérie Orsat. "Microwave- Assisted Extraction of Flavonoids: A Review", Food and Bioprocess Technology, 2011 Publication	<1%
13	www.soft-tox.org	<1%
14	www.tandfonline.com	<1%

Queiroz, Emerson F., Kurt Hostettmann, and Jean-Luc Wolfender. "Modern Approaches in the Search for New Active Compounds from Crude Extracts of Natural Sources", Plant Bioactives and Drug Discovery Principles Practice and Perspectives, 2012. Publication

<1%

16 Jean-Luc Wolfender, Jean-Marc Nuzillard, Justin J. J. van der Hooft, Jean-Hugues Renault, Samuel Bertrand. "Accelerating Metabolite Identification in Natural Product Research: Toward an Ideal Combination of Liquid Chromatography–High-Resolution Tandem Mass Spectrometry and NMR Profiling, Databases, and Chemometrics ", Analytical Chemistry, 2018 Publication

17	www.omicsonline.org	<1%
18	Submitted to University of Huddersfield Student Paper	<1%
19	www.medicinainternacba.com.ar	<1%
20	repository.unhas.ac.id	<1%

documents.mx

21

22	Submitted to iGroup Student Paper	<1%
23	Green Chemistry and Sustainable Technology, 2014. Publication	<1 %
24	Submitted to University of Strathclyde Student Paper	<1%
25	 FA Weltzien, G.I Hemre, J.O Evjemo, Y Olsen, H.J Fyhn. "β-Hydroxybutyrate in developing nauplii of brine shrimp (Artemia franciscana K.) under feeding and non-feeding conditions", Comparative Biochemistry and Physiology Part B: Biochemistry and Molecular Biology, 2000 Publication 	<1%
26	G. Brusotti, I. Cesari, A. Dentamaro, G. Caccialanza, G. Massolini. "Isolation and characterization of bioactive compounds from plant resources: The role of analysis in the ethnopharmacological approach", Journal of Pharmaceutical and Biomedical Analysis, 2014 Publication	<1%

Exclude quotes	On	Exclude matches	Off
Exclude bibliography	On		

Metabolite Profiling Analysis of Conventional and Non-Conventional Extraction Methods on Secondary Metabolite from Peperomia pellucida (L.) Kunth using UPLC-QToF-MS/MS System

GRADEMARK REPORT	
FINAL GRADE	GENERAL COMMENTS
PAGE 1	
PAGE 2	
PAGE 3	
PAGE 4	
PAGE 5	